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DW DENTAL
WORLD
ITALIAN MANUFACTURER

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


An Initial Investigation of the Bending and Torsional Properties of Nitinol Root Canal Files

Harmeet Walia, BDS, MDS, MS, MS, William A. Brantley, BS, MS, PhD, and Harold Gerstein, BS, DDS

Root canal files in size #15 and triangular cross-sections were fabricated from 0.020-inch diameter arch wires of Nitinol, a nickel-titanium orthodontic alloy with a very low modulus of elasticity. A unique manufacturing process was used in which the fluted structure of a K-type file was machined directly on the starting wire blanks. The Nitinol files were found to have two to three times more elastic flexibility in bending and torsion, as well as superior resistance to torsional fracture, compared with size #15 stainless steel files manufactured by the same process. The fracture surfaces for clockwise and counterclockwise torsion were observed with the scanning electron microscope and exhibited a largely flat morphology for files of both alloy types and torsional pretesting modes. It was possible to permanently precurve the Nitinol files in the manner often used by clinicians with stainless steel files. These results suggest that the Nitinol files may be promising for the instrumentation of curved canals, and evaluations of mechanical properties and *in vitro* cutting efficiency are in progress for size #35 instruments.

It is well known by clinicians that inadvertent procedural errors can occasionally arise during the instrumentation of curved canals. These misfortunes include ledge or zip formation, perforation of the canal, and separation or fracture of the instrument (1). As a consequence, the root canal morphology is adversely altered, a violation of the basic principle that endodontic preparation is to retain the original shape of the canal. Clinicians have adopted various methods to circumvent problems with the preparation of curved canals, such as precurving instruments and using a telescopic filing technique (1-3). Weine (4) has suggested that clinicians might remove the tips of instruments at chairside to make intermediate sizes for use in the preparation of curved canals.

The procedural errors which may occur during the instrumentation of curved canals have a common genesis: the basic stiffness of the stainless steel alloys (5) utilized for the manufacture of root canal files and reamers. Moreover, there is a substantial rise in instrument stiffness with increasing instrument size (6). For example, with the stainless steel files and reamers, the smaller sizes of instruments have considerably

greater flexibility and can conform much better to the morphology of curved canals.

While manufacturers have recently marketed a number of new instruments based upon different cross-sectional shapes, design concepts, and fabrication procedures, in a quest for improved cutting efficiency (7) and flexibility (8), all of these brands have been fabricated from stainless steel. In this article we report the first use of an entirely new metallurgical system, the Nitinol nickel-titanium orthodontic wire alloy (9), for the fabrication of endodontic files. The Nitinol alloy has a very low modulus of elasticity, only one-fourth to one-fifth the value for stainless steel, and a very wide range for elastic deformation.

The purposes of this initial study were to investigate the feasibility of manufacturing root canal files from Nitinol and to evaluate the bending and torsional properties of these instruments. The results of our laboratory study suggest the possibility of a new generation of files, possessing a degree of flexibility which may be ideally suited for instrumenting curved canals.

MATERIALS AND METHODS

Standard preformed Nitinol arch wire blanks, 0.020 inch in diameter, were obtained (Unitek Corp., Monrovia, CA), and two 2-inch straight segments from each arch wire were used for instrument fabrication. A unique file manufacturing process was used (Quality Dental Products, Johnson City, TN), in which the fluted cross-sectional shape was machined directly on the wire blank, rather than the conventional (10) manufacturing procedure of twisting the ground and tapered blank. For this initial feasibility study, experimental Nitinol root canal files were fabricated in size #15 and triangular cross-sections, for comparison to size #15 stainless steel files with the same cross-sectional shape and manufactured by the same process, which served as the controls.

The Nitinol and stainless steel files were evaluated in three mechanical testing modes of cantilever bending, clockwise torsion, and counterclockwise torsion, following the experimental methods previously used by Krupp et al. (8). Values of bending and torsional moment were measured with a sensitive torque meter (model 783-C-1; Power Instruments, Inc., Skokie, IL), using a manual-loading experimental procedure and an apparatus based upon the original form of American Dental Association specification no. 28 (11). All specimens were subjected to bending or twisting at a point 3

The forms of the bending curves in Fig. 5 indicate that permanent deformation of the 3-mm apical regions of the stainless steel files began at a bend angle of approximately 30 degrees, but that the apical regions of the Nitinol files were undergoing largely elastic deformation even at bend angles of 90 degrees. The latter was supported by visual observations of the Nitinol files after unloading, where very little, if any, permanent bends were evident.

The Nitinol files also exhibited considerably greater resistance to fracture in torsion than the stainless steel files. For


Fig. 6. Clockwise torsion test results for the size #15 Nitinol and stainless steel files.


Fig. 7. Counterclockwise torsion test results for the size #15 Nitinol and stainless steel files. The two initial data points for the Nitinol files were determined with the torque meter, and the two plots were drawn to intersect the origin. Both of these considerations are pronounced in Figs. 5 and 6.

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
MOST OF the NEW FILES
IN THE MARKET
ARE HEAT TREATED

*Adapting The Anatomy
to the File Sequence ?*

Or

*Adapting The File Sequence
to the Anatomy ?*

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
alkal
canal sea

konstantinos Sidiropoulos , Elisabeth Koulaouzidou ,
Economides 

of Dental Tissues, School of Health Sciences, Faculty of Dentistry,
Athens, Greece

properties (pH and
hydroxi

coloration
presence of sodium

Angélica Marciano¹ · Marco Antonio Hu

Received: 11 November 2014 / Accepted: 18 March 2015
Springer-Verlag Berlin Heidelberg 2015

The aim of this research was to analyze the coronal aspect of the root canal *ex vivo* and the effect of red blood cells on this discoloration. Canals were prepared from the apical aspect with either wMTA + saline ($n = 18$), or controls ($n = 4 + 4$) (blood + saline). Discoloration was assessed according to the CIE L*a*b* color system using standardized digital photographs at baseline, day 1, and day 7. The results were analyzed by using a two-way ANOVA test. The results showed that the presence of sodium bismuth oxide in the wMTA aggregate caused a significant increase in discoloration over time.

ology

Root Canal Tooth Discoloration by a Silicate-based Material

Daniel Felman, BDS_c, DCD, and Peter Paraschos, BDS_c, DCD

Abstract
Introduction: This study assessed and characterized the coronal aspect of the root canal *ex vivo* and the effect of red blood cells on this discoloration. Canals were prepared from the apical aspect with either wMTA + saline ($n = 18$), or controls ($n = 4 + 4$) (blood + saline). Discoloration was assessed according to the CIE L*a*b* color system using standardized digital photographs at baseline, day 1, and day 7. The results were analyzed by using a two-way ANOVA test. The results showed that the presence of sodium bismuth oxide in the wMTA aggregate caused a significant increase in discoloration over time.

Light and
Silicate-based Material

allés, DDS, MsC,* Montse Mercadé, DDS,
Bourdelande, BSc, PhD,[†] and Miguel Roig

Abstract
Introduction: Difficult handling, long setting time, and discoloration are important drawbacks of white mineral trioxide aggregate (wMTA). The development of a new calcium silicate-based material (CSM), has overcome some of these drawbacks. However, there are no available studies on the effect of light irradiation in an *in vitro* model. The present study evaluated the effect of light irradiation on the discoloration of wMTA and CSM. The results showed that the presence of sodium bismuth oxide in the wMTA aggregate caused a significant increase in discoloration over time.

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OUTCOME

Outcome of Non-Surgical Root Canal Treatment Using a Single-cone Technique with a Sequence Bioceramic Sealer: A Retrospective Analysis

Elizabeth A. Chybowski, DDS,* Gerald N. Glickman, DDS, MS,† Eric Fleury, DDS, MS,‡ Eric Solomon, DDS, MS,‡ and Jianing J...

Abstract
Introduction: One of the important steps in root canal treatment is to create a well-sealed root canal system. Sequence BC Sealer (BC; Brasseler USA, Savannah, GA) has several beneficial properties and thus has been incorporated into the practitioner's armamentarium. No study to date have evaluated the clinical success of using a single-cone technique with a minimum of treatment factors were...

Conclusion: An important goal of root canal treatment is to properly seal the canal system after cleaning and shaping. However, irregularities such as fins, isthmuses, and lateral canals are often present and can pose challenges to clinicians during obturation. Anatomical spaces can have treatment (5). Historically, root canal obturation with poor root canal obturation techniques were...

3-month follow-up of primary root canal treatment of teeth obturated with a hydraulic sealer

Giulia Bardini¹ • Laura Casula² • Emanuele Ambu¹ • D...

Received: 3 February 2020 / Accepted: 15 September 2020 / Published online: 15 September 2020
© The Author(s) 2021, corrected publication 2021

Abstract
Objectives This randomized, controlled, pilot study compared the outcomes of root canal obturation either with a novel bioactive sealer and warm vertical compaction or with a hydraulic sealer and warm vertical compaction.

Materials and methods Sixty-nine patients with root canal treatment were included in the study. The technique with BioRoot™ RCS (Septodont, Brétigny-sur-Orge, France) was compared with the hydraulic sealer technique with BioRoot™ RCS (Septodont, Brétigny-sur-Orge, France) in a randomized controlled trial.

Outcome of Root Canal Treatment with a Calcium Silicate Root Canal Sealer: A Non-Randomized Clinical Study

Zavattini^{1,*}, Alan Knight¹, Federico Foschi¹, ...

Abstract
Objectives This study aimed to evaluate the clinical outcome of root canal treatment with a calcium silicate root canal sealer compared with a bioceramic-based sealer.

Materials and methods Sixty-nine patients with root canal treatment were included in the study. The technique with a calcium silicate root canal sealer was compared with the technique with a bioceramic-based sealer in a non-randomized clinical study.

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